**Technical Specification**

Part No.	Head Alloy	Product Weight	Max Patient Weight	Torque
3CT5	Titanium	26gr	80kg/175lbs	M6x10 – 12Nm / M4x10 – 9Nm
3TC1	Titanium	44gr	80kg/175lbs	M6x10 – 12Nm / M4x10 – 9Nm
3TC2	Titanium	51gr	80kg/175lbs	M6x10 – 12Nm / M4x10 – 9Nm

**Description and Purpose**

These instructions are for use by the practitioner. This product is designed to be used for lower limb external prosthetic length, toe direction and alignment corrections. Intended for single user only. This product can be used for mobility classes 1-4. There are no known contraindications if used in accordance with these instructions. Ensure that the user has understood any instructions for use, drawing particular attention to the maintenance information.

**Construction**

Use 3mm Allen wrench for all screws. Before fitting this product onto the patient, all parts must be joined properly and any set screws must be secured enough for the patient to walk. In order to make final adjustments, there should be no loose screws. After the final adjustments, apply single drop of removable Loctite® 242 to each of the set/cap screws and then torque the product.

**Function**

Connects directly to socket (4TC7 lamination plate) with M5 countersunk screws (10Nm, not provided). Use 3mm Allen wrench for set screws. No Loctite recommended for the 5mm countersunk mounting screws. Rotation set M4x10 screw torque 9 Nm (3TC1, 3TC2) M6x10 angular adjustment set screw torque 12 Nm. Improper installation or uneven, insufficient or excessive torque may cause bolt or adapter failure. Always use quality metric tools for installation. When tightened securely using Loctite on the pyramid screws only and torque settings shown, the device offers a secure way to connect two parts of the limb. Once the desired toe direction is achieved, torque the screws to 9Nm. Use one drop of removable Loctite 242 after the final adjustment. Improper installation or uneven, insufficient or excessive torque may cause bolt or adapter failure. Always use quality metric tools for installation. When tightened securely using Loctite on the pyramid screws only and torque settings shown, the device offers a secure way to connect two parts of the limb.

**Limitations on Use**

**Intended life:** A local risk assessment should be carried out based upon activity and usage.  
**Lifting loads:** User weight and activity is governed by the stated limits. Load carrying by the user should be based on a local risk assessment.

**Environment:** Thoroughly rinse with fresh water after use in abrasive environments such as those that may contain sand or grit, for example, to prevent wear or damage to moving parts. Thoroughly rinse with fresh water after use in salt or chlorinated water. Exclusively for use between -15°C and 50°C (5°F to 122°F).

**Maintenance**

Maintenance is to be carried out by competent personnel only. An annual visual inspection is recommended. During maintenance, it is appropriate to consider any weight or activity increase by the user. If the prosthesis is exposed to circumstances or impact for which it was not designed, inspection by a qualified practitioner is recommended.

- Ensure all screws are secure. If not, remove and clean screws, reapply Loctite to the pyramid screws only, tighten all screws to the correct torque settings.
- Check for visual defects that may affect proper function.
- The user should be advised to contact their practitioner if:
  - Their condition changes.
  - Performance of the device changes.
  - Changes in performance may include:
    - o Any unusual noises
    - o Excessive play or loss of alignment

**Warranty**

All MEDEX Modular Adapters are unconditionally guaranteed for 2 years. Cutting, drilling, modifying or using the component other than as described in this manual will void the warranty.

**Liability**

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device.

**CE Conformity**

This product meets the requirements of 93/42/EEC guidelines for medical products. This product has been classified as a Class 1 Product according to the classification criteria outlined in Appendix IX of the guidelines.

MEDEX INTERNATIONAL, INC. – 2720 Duvall Road, Burtonsville, MD 20866 – Phone: 1-888-886-2420 – Fax: 301-657-2796 – Last Revised 4/1/25

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Part No.	Material	Gewicht	Gewichtslimit	Drehmoment
3CT5	Titan	26gr	80kg/175lbs	M6x10 – 12Nm / M4x10 – 9Nm
3TC1	Titan	44gr	80kg/175lbs	M6x10 – 12Nm / M4x10 – 9Nm
3TC2	Titan	51gr	80kg/175lbs	M6x10 – 12Nm / M4x10 – 9Nm

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