

## Waterproof Structural Components



**Components include grade A2 stainless steel corrosion resistant screws**

### **Warranty:**

Medex components below are guaranteed to be waterproof in chlorinated, salt, or fresh water environments for up to 12 months, plus an additional 12 months if rinsed with fresh water after each use in chlorinated or salt water.

Warranty covers all grade A2 stainless steel screws.



Part Name	Part #
SACH FOOT ADAPTER [8MM BOLT]	2T09
SACH FOOT ADAPTER [10MM BOLT]	2T10
SACH FOOT ADAPTER [ROUND, 10MM BOLT]	2T11
LOW PROFILE PYLON WITH RECEIVER - 400 x 30 x 2.35MM	2T12
HEAVY DUTY TUBE ADAPTER [RATED 450 LBS]	2T13
LONG TUBE WITH PYRAMID ADAPTER - 400 x 30 x 2.5MM	2T14
LONG TUBE WITH 4 SCREW ADAPTER - 400 x 30 x 2.5MM	2T15
SHORT TUBE WITH 4 SCREW ADAPTER - 250 x 30 x 2.5MM	2T16
EXTRA LONG TUBE W/ 4 SCREW ADAPTER - 450 x 30 x 2.35MM	2T17
TUBE CLAMP ADAPTER	2T20
TUBE CLAMP ADAPTER [10 DEGREE ANGLE]	2T22
HEAVY DUTY TUBE CLAMP ADAPTER [RATED 450 LBS]	2T23
PYRAMID ADAPTER W/ 4 HOLES [10MM CENTER HOLE]	3T09
PYRAMID ADAPTER W/ 4 HOLES	3T10
PYRAMID ADAPTER W/ 4 HOLES [8.2MM CENTER HOLE]	3T10
PYRAMID ADAPTER W/ 4 HOLES [AXIAL OFFSET]	3T11
PYRAMID ADAPTER W/ 4 HOLES [DIAGONAL OFFSET]	3T12
PYRAMID ADAPTER W/ 4 PRONGS [HEAVY DUTY]	3T13
PYRAMID ADAPTER W/ 4 PRONGS	3T14
PYRAMID ADAPTER W/ 4 PRONGS - ROTATABLE	3T15
PYRAMID RECEIVER W/ 4 PRONGS - ROTATABLE	3T16
PYRAMID RECEIVER W/ 4 PRONGS	3T17
HEAVY DUTY PYRAMID ADAP. W/4 HOLES 8.2MM CENTER HOLE [THREADED ROTATION]	3T18
HEAVY DUTY PYRAMID ADAP. W/4 HOLES 10MM CENTER HOLE [THREADED ROTATION]	3T19
4 SCREW PYRAMID RECEIVER	3T20
HEAVY DUTY PYRAMID RECEIVER W/ 4 HOLES [THREADED ROTATION]	3T23
PYRAMID ADAPTER W/ 4 HOLES [THREADED ROTATION][DIAGONAL OFFSET]	3T24
PYRAMID ADAPTER W/ 4 HOLES [THREADED ROTATION][AXIAL OFFSET]	3T25
PYRAMID RECEIVER W/ 4 HOLES [THREADED ROTATION] [DIAGONAL OFFSET]	3T26
PYRAMID RECEIVER W/ 4 HOLES [THREADED ROTATION] [AXIAL OFFSET]	3T27
8 SCREW DOUBLE HEAD ADAPTER - 32MM	3T50
8 SCREW DOUBLE HEAD ADAPTER - 45MM	3T51
8 SCREW DOUBLE HEAD ADAPTER - 60MM	3T52
8 SCREW DOUBLE HEAD ADAPTER - 75MM	3T53
ROTATABLE DOUBLE FEMALE STAT ADAPTER 34MM	3T54
ROTATABLE DOUBLE FEMALE STAT ADAPTER 42MM	3T55
ROTATABLE DOUBLE FEMALE STAT ADAPTER 53MM	3T56
ROTATABLE DOUBLE FEMALE STAT ADAPTER 66MM	3T57
4 SCREW AND PYRAMID DOUBLE HEAD ADAPTER	3T60
ROTATABLE MALE/FEMALE STAT ADAPTER 42MM	3T61
ROTATABLE MALE/FEMALE STAT ADAPTER 34MM	3T63
HEIGHT ADJUSTABLE & ROTATABLE STAT ADAPTER 50MM to 57MM	3T65
HEIGHT ADJUSTABLE & ROTATABLE STAT ADAPTER 58MM to 64MM	3T66
BOTH END PYRAMID ADAPTER [SHORT]	3T70
TITANIUM 4 HOLE ROTATOR BASE	3T90
RECEIVER INSERT	4T08
PYRAMID INSERT	4T09
PYRAMID INSERT 2MM LONGER	4T10
PYRAMID INSERT 4MM LONGER	4T11
3-ARM LAMINATION ANCHOR FRAME	4S35
3 ARM A/K ROTATABLE LAMINATION ANCHOR [W/ TITANIUM RECEIVER]	4ST6
3 ARM A/K ROTATABLE LAMINATION ANCHOR [W/ TITANIUM PYRAMID]	4ST8

---

Quality Products for the Prosthetic Industry

MEDEX INTERNATIONAL, INC. - PO Box 21 Kensington, MD 20814 – Phone: 1 888-886-2420 – Fax: 301-657-2796